

<b>pts panels<sup>®</sup></b> Multi-Chemistry Controls    Range Card	Level 1		Level 2	
	mg/dL	mmol/L	mg/dL	mmol/L
<b>Lot #:</b> <b>MC41</b> <b>Expiration Date:</b> <b>2023-05-31</b>				
<b>Single-Analyte Test Strips</b>				
Cholesterol	130-230	3.37-5.96	140-305	3.63-7.90
HDL Cholesterol* <i>(CardioChek<sup>®</sup> Plus analyzer version 1.09 or above)</i>	≤25-60	≤0.65-1.55	25-75	0.65-1.94
Triglycerides	90-200	1.02-2.26	160-305	1.81-3.45
Glucose/eGLU	40-120	2.22-6.66	170-450	9.44-24.98
<b>Multi-Analyte Test Strips</b>				
Cholesterol	105-245	2.72-6.35	150-350	3.89-9.07
HDL Cholesterol* <i>(CardioChek<sup>®</sup> Plus analyzer version 1.09 or above)</i>	≤15-60	≤0.39-1.55	25-85	0.65-2.20
Triglycerides	80-205	0.90-2.32	150-310	1.70-3.50
Glucose	25-130	1.39-7.22	130-410	7.22-22.76
* For CardioChek Plus (version 1.08 or below), CardioChek PA, and CardioChek analyzers, HDL Cholesterol requires separate HDL Cholesterol Controls (REF 722). PTS Diagnostics Customer Service: 877-870-5610   Fax: 317-870-5608   Outside of US: +1-317-870-5610 Email: customerservice@ptsdiagnostics.com				
				PS-002519 Rev. 14